
MATR X® **MAC BACK**

ADJUSTABLE CONTOUR









Introducing a whole new approach to bespoke back support, the Matrx MAC Adjustable Contour Back Series! Completely configurable for the widest range of body shapes and conditions, and adjustable for changes,

the MAC back offers a whole new way to provide custom fit, positioning and comfort!















MAC - HUG MODEL



Soft, conformable StaminaFibre® for maximum immersion and envelopment



- StaminaFibre® insert, machine wash and dry friendly
- Highly immersive micro climate control
- Recommended for accommodation of mild postural asymmetries



 Customisable shape interface for accommodation of orthopedic asymmetries



MAC - HUG MODEL/SHAPE



SHAPE Insert layered under HUG StaminaFibre®



- Positioning interface with configurable HR foam pieces
- Easily adjustable for future changes
- Recommended for accommodation of moderate postural asymmetries

MAC - DUAL LAYER FOAM MODEL





Supersoft/firm HR foam layers for comfort and support






- Supersoft over firm HR double layer foam
- Support and comfort even in tilt/recline applications
- Recommended for symmetrical postural presentation







BACK SIZES

WIDTH

HEIGHT		16"	18"	20"
	FITS W/C	15"-18"	17"-20"	19"-22"
	16"		•	
	18"	•	•	•
	20''	•	•	•
	22"		•	

· Up to 133kg weight capacity

OPTIONAL HARDWARE



- 1. Extended Range EasySet
- Provides 4" forward or rear depth adjustment
- 2. Fixed EasySet
 Recommended for tilt/recline applications
 Requires tools for removal
- 3. Heavy-Duty EasySet
- Recommended for users with high tone or heavy intended use
 Split
- Allows clearance of frame obstacles, quick release



PostureFlex

Comfortable, quick release anterior trunk (abdominal) supports. Available in 2 sizes.



BackBone

Contoured PSIS/Sacrolumbar Support.



Quad Release Tab

Quad release tab for PostureFlex.



Matrx Privacy Flaps

Easy attachment to all Matrx Backs with existing sewn on hook & loop. Available in 4 sizes.



PostureFlex Tether Strap

Required to mount PostureFlex. Compatible with MAC - 3.5" Contour Laterals.



Chest Harness Interface

Allows precise placement of Shoulder Strap Guides for Chest harnesses. Available in 2 sizes.





Canada 84 Citation Dr., Unit I Concord, ON L4K 3CI (866) 748-7943 (905)-695-0 | 34 Fax: 905-695-0138

© 2024 Matrx is a registered trademark of Motion Concepts. Specifications are subject to change without notification.

TRD0733 Rev A

0800212333 or info@c1south.co.nz 1/7 Newton Place, Hamilton 3204