

PFLEX-MAC, PostureFlex Abdominal Support FOR MATRX MAC BACKS

1

6.5" (165 mm)

PFLEX-MAC-M
Adjustment Range:
16 - 24" (41 - 61cm)

2

7" (178 mm)

PFLEX-MAC-L
Adjustment Range:
20 - 34" (51 - 86cm)

3

4

5

6

(Binding Nut faces outward)


Item	Qty	Description
1	1	PFlex-M Abdominal Support- 6.5 H (165mm)
2	1	PFlex-L Abdominal Support- 7 H (178mm)
3	2	10-24 x 3/8" Standard Truss Machine Screw
4	2	3/16" Flat Washer, SAE, BZ (M5 Compatible)
5	2	10-24 x 1/4" Binding Nut
6	2	Tether Strap w/ Slider and Mounting Clip

Loop

Slider

Tether Strap Routing



TETHER STRAP INSTALLATION:






*Note: The Plastic Cap can be re-installed following tether strap installation/adjustment

PFLEX-MAC INSTALLATION:

REAR VIEW



FRONT VIEW



Motion Concepts

505 Cityview Blvd, Unit 1
Vaughan, Ontario, L4H 0L8, Canada
Tol Free: (1) 866.748.7943
Tel: (1) 905.695.0134

Motion Concepts USA:

700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
Tol Free: (1) 888.433.6818